

**HACH's answers to
wastewater challenges
PART VIII : Self
monitoring**



Disclaimer

Hach webinars are educational in nature and not an endorsement of other manufacturers' treatment technologies or processes.

This document contains information owned by Hach. Personal use, copying or distribution is allowed for webinar participants for their own internal purposes, when mentioning Hach as the source.



Speakers



Carsten SCHULZ

Product Applications Manager




Philippe PONS

Application Development Manager

WWTP monitoring

Overview of Hach solutions





Agenda

1. Laboratory measurements: Why?
2. Automatic water samplers
3. Laboratory: Instruments and analytical techniques
 - Spectrophotometry
 - Benefits of cuvette tests
 - COD
 - TOC
 - Phosphate
 - Nitrogen parameters



Laboratory measurements: Why?

Main goal of Municipal plants is to correctly treat water to protect environment and meet regulatory requirements

- On-line instruments are used to monitor and optimize the efficiency of the process in real time
- Laboratory instruments are used to measure punctually samples:
 - To monitor the process if there is not online instruments
 - To calibrate / check the online instruments
 - For self-monitoring in accordance with EU / Local regulation

Only laboratory results are accepted by authorities to prove the proper functioning of the plant (self-monitoring)!

Self-monitoring process for Muni WW

Urban Wastewater Treatment Directive 91/271/EEC (UWWTD) - 21/05/91

(Applicable for plants > 2000 PE)

Defines...

- The parameters to be monitored
DBO5 / COD / TSS / TP / TN
- The analytical technique (not mandatory)

... And, depending on the plant size and on the aera type (sensible or not),...

- The maximum limit for each parameters
- The frequency of self-measurements

Local regulation can be stricter than EU recommendations

**The UWWTD is under revision!
New version will be more restrictive**

Measurement frequency

Minimum number of samples to be taken at regular intervals over a full year and according to the size of the treatment plant:

- P.E. between 2,000 and 9,999:
 - 12 samples in the first year.
 - 4 samples in subsequent years if it can be demonstrated that the water complies with the provisions of this Directive in the first year; if any of the 4 samples do not meet the standards, 12 samples are taken the following year.
- P.E. between 10 000 and 49 999: 12 samples / year.
- P.E. of 50 000 or more: 24 samples / year.

(PE: Organic load of DBO5 = 60 grs/day)



Laboratory analysis: 5 Steps

- Sampling
- Sample Preservation
- Sample Preparation
- Use of Standards (check and calibration)
- Analysis



The analysis is only as good as the sample!




Automatic Water Samplers



SAMPLING METHODS: How and when?

HOW?: ISO Standard: ISO 5667-10 :2020

WHEN?: Local regulation



Ex: Fench water agencies

Critères de représentativité <i>et valeur à respecter</i>	Fréquence de vérification préconisée en fonction de la fréquence d'utilisation du préleveur
<u>Vérification lors d'un bilan 24 h :</u> <ul style="list-style-type: none"> • Volume unitaire d'un prélèvement (1) <i>≥ 50 ml</i> • Répétabilité ou distribution des volumes (2) <i>Ecart maximum toléré ≤ 5%</i> • Nombre de prélèvement en 24h <i>≥ 100</i> • Volume total prélevé en 24h <i>≥ 5 litres</i> 	Mensuelle	1 fois par jour
	Trimestrielle	Entre 1 fois par jour et 1 fois par mois
	Semestrielle	Moins d'1 fois par mois
<u>Vérification lors d'un bilan 24 h :</u> <ul style="list-style-type: none"> • Température de l'enceinte <i>3°C ± 2°C (NF ISO 5667-3)</i> ou <i>2°C ± 2°C sans congélation (NF ISO 5667-10)</i> ou <i>4°C ± 2°C (règles de l'art)</i> 	Au moins 2 fois par an	Toutes fréquences confondues
<ul style="list-style-type: none"> • Vitesse d'aspiration <i>≥ 0,5 m/s</i> 		



SAMPLING METHODS

- **Manual Sampling**



- **Automatic Sampling**



Manual sample: SPOT CHECK




Quantity
 Q



Sample
Volume
per
Sample



Sample
Volume
per
Bottle



AUTOMATIC WATER SAMPLING




- **Sampling methods**


- Time related
- Connected to flow meter
 - Constant Volume/Variable Time (CV/VT)
 - Constant Time/Variable Volume (CT/VV)

- **2 sampling techniques**


- Peristaltic pumps
- Pressure vacuum



TIME RELATED SAMPLE



Automatic water sampler with connected flow meter




Sampler takes in consideration the sample flow

2 sampling methods


- Constant Volume/Variable Time (CV/VT)
- Constant Time/Variable Volume(CT/VV)



FLOW PROPORTIONAL SAMPLE: Constant Volume/Variable Time (CV/VT)




FLOW PROPORTIONAL SAMPLE: Constant Time/Variable Volume(CT/VV)



Why to cool samples?


Cooling samples at 4°C inhibits bacteria activity → sample content won't be modified




Active or passive cooling available

Technologies, how to take samples automatically?

Peristaltic technology (Sigma)



Pressure vacuum (Buehler)



What says the norm regarding Peristaltic pump?



ISO 5667-10:2020
Annex F

Table F.1 — Advantages and disadvantages of two main types of pumping

Pumping type	Advantages	Disadvantages
Peristaltic pump (PP)	Only the desired volume of sample is pumped	Accuracy of unit volume: difficult to obtain especially for small volumes (<100 ml)
	Possible rinsing with sample to avoid cross contamination	Checking the unit volume required over time
	Sample contact only with the inside of the suction pipe and the liquid detector	Possibility of deposits in pipes which can lead to cross contamination
	No sedimentation of solids because sample is in permanent movement from the sampling point to the collecting flask	Unstained solid integrity due to pipe crushing (less representative sample for suspended solids, possible biological modification)
	Cleaning and preventive maintenance low, easy and less expensive	Internal diameter limited to 9,5 mm
	Crushing pipe and rollers are the only wear parts. Easy to change and low cost	Suction hard solids (sand) can damage the pump and result in high maintenance costs
	Better overall behaviour under particular sampling conditions (hot fluids, fouling, presence of suspended matter)	Frequent exchange of the crushing pipe. Risk of rupture if rarely replaced and / or in the presence of hard solids in the sample.
	Slower fouling of elements in contact with the fluid	
	High suction heights	



Peristaltic pump samplers = HACH SIGMA

**All Weather
Refrigerated (AWR)**



Refrigerated

Portable



What says the norm regarding Vacuum pump



ISO 5667-10:2020
Annex F



Table F.1 (continued)

Pumping type	Advantages	Disadvantages
Vacuum pump (VAP)	Repeatable and accurate sample volumes	Possible sedimentation and degassing in the bowl during adjustment to the requested volume
	Easy calibration of the unit volume	Risk of clogging the pressure sensor
	Reduced risk of cross-contamination by pre- and post-sampling purges	Mandatory detection of the arrival of the sample
	Solids integrity	Risk of permanent deposition in the bowl (e.g. grease) which can alter the volume and composition of the sample
	Pipe diameter may be greater than 9 mm	Various materials in contact with the sample
	Easy cleaning of the bowl (disassembly)	More complex system
	Less frequent change of wear parts (vacuum pump diaphragm, pinch valve pipe)	Cleaning or rinsing of pipes more difficult to carry out "in-line" without disassembly (problem of absence of conductivity when rinsing with ultra-pure water)
	Less Energy Consumption	
Possible rinsing with sample to avoid cross contamination	Limited pumping height (7 m)	



Pressure/Vacuum samplers = BÜHLER



How to select the correct sampler?

STATIONARY

✓ Conditions of sampling

Is the sampling done on loaded pipe or at atmospheric pressure?

✓ Numbers of bottles

Which bottle do you need?

✓ Bottle material

Do you prefer plastic or glass bottles?

✓ Refrigerated model (portable)

Do you need refrigerated model?

✓ STAINLESS STEEL ou PE

What materials do you want to take it?

✓ Presence of H₂S

Is there H₂S?

If yes, choose reinforced stainless steel 316/GF




1 x 25L



4 x 14L



24 x 1L



How to select the correct sampler?

PORTABLE

- ✓ Peristaltic pump or vacuum technique

Is the sampling done on loaded pipe or at atmospheric pressure? At which height?

- ✓ Numbers of bottles



Which bottle do you need?

- ✓ Bottle material

Do you want plastic or glass bottles?

- ✓ Insulated or refrigerated model

Do you want an insulated or refrigerated model?



Laboratory Instruments and analytical technique



Analytical Technologies

Photometric (Light)

- Visual
- **Spectrophotometry (Colorimetry)**
- Nephelometry (Turbidity)
- Luminescent (LDO)

Electrochemical (Probe)

- Potentiometric (pH, ISE...)
- Polarographic (DO)
- Conductometric
- ORP

Titrametric

- Drop Count
- Digital Titrator
- Traditional Buret
- Electrochemical titration (Probes)

Microbiological

- Qualitative
- Quantitative


Gravimetric

- Total Dissolved Solids
- Total Suspended Solids



Spectrophotometry



Principle of photometry



- Water parameters do not have an apparent color
- When adding specific reagents, specific colors develop



- The colored sample absorbs light
- The absorption of the light is proportional to the concentration of the parameter

Photometric System

Designed as a system

Hach chemistries and instruments are designed to work as an integrated system, eliminating calibration and variability errors from the outset so you can be certain of your results.

Reagent stability

Formulations and packaging are designed to provide long product life to ensure accurate results after shipping and storage and reduce the effects of common sample interferences.

Support

We offer a dedicated support team who are trained to solve your technical problems, backed by decades of experience.





Traditional Hach Tests

- Bulk liquids
- AccuVac Ampules
- Powder Pillows
- PermaChem powder reagents
- **LCK cuvette tests**
- Chemkey



Why Perform Self Monitoring with Cuvette



Benefits of cuvette tests

- Safety
 - No contact with chemistries
 - Less chemical reagent than other methods
- Easy to use & ideal for wastewater parameters
- Rapid and accurate results
- Easy updates
- Environmental protection (recycling)
- Standardized Methods



Safety

No contact with chemistries

- Cap with reagent
 - DosiCap / DosiCap Zip for freeze-dried reagents
 - No powder or liquid reagent handling with several tests
- Cuvette pre-filled reagent
- No waste handling
- No glassware cleaning




Safety

~95% fewer chemical reagents than traditional standard method

COD with
Standard Method



Hach COD Test



Easy to use

- Complete solution & ready to use
- All required reagents in one box
- Ideal for wastewater parameters (digestion)
- Easy-to-understand and structured working procedure
 - Pictograms inside box lid
 - Full procedure available by download
- Quality certificate by download

Quality certificate
 Technical data for Validation
 of LCK400 (0-1000 mg/l COD
 in conformity with ISO 15705)


Quality certificate

Technical data for cuvette test LCK400

Sensitivity	0.0004 Abs./(mg/l)
Ordinate intersect	0.003 Abs.
Residual standard deviation	0.0007 Abs.
Method variation coefficient	0.32 %
Method standard deviation	1.7 mg/l
Confidence interval (95%)	± 4.2 mg/l

Technical data in conformity with DIN 32645

Detection limit	6.1 mg/l
Quantitation limit	21.0 mg/l



The technical data for cuvette test LCK400 were determined in conformity with ISO 8466-1 and DIN 38402 A51 „Calibration of analysis methods“.

The series of the smallest and largest calibration standards exhibit normal distribution and are outlier- and trend-free. The calibration gives a linear function.

The detection and the quantitation limits were determined in conformity with DIN 32645.


Result	Confidence interval
200 mg/l	± 4.5 mg/l
400 mg/l	± 4.3 mg/l
600 mg/l	± 4.2 mg/l
800 mg/l	± 4.3 mg/l
1000 mg/l	± 4.6 mg/l

HACH LANGE GmbH
 Quality Management


R. Kloos
 Dr. Ralf Kloos




Procedure




1. Carefully pipet 1.0 mL of sample.



2. Carefully pipet 0.2 mL of solution A.




3. Close the cuvette and invert a few times until no more streaks can be seen.



4. After 15 minutes, thoroughly clean the outside of the cuvette and evaluate.




Accurate Results – 10fold measurement




- **Reference beam** technology corrects net deviation or lamp aging effects
- **10-fold-measurement** (including averaging and outlier elimination) to avoid wrong results from scratches, dirt and fingerprints
- **Warning** in clear text messages when turbidity or dirty cuvettes in place

Accurate Results – 2D code



What comes with the 2D-code:

- **Test #** and automated check if right test data are available
- **Test Expiration date:** Message if reagents are expired
- **Lot number** (traceability)
- **Lot specific factor** (Truecal)



Accurate Results – Truecal


The problem: Inaccurate & Unreliable Testing

- 8-10% of labs in the EU fail proficiency testing on an annual basis
- Bad readings and potential limit violations are not identified
- Failing proficiency testing with potential loss of laboratory accreditation

The solution:

Truecal™

- Lot after lot - Be Right
- Truecal with each vial includes the calibration data for each individual lot, reducing variation in results.
- This allows to meet reporting results and to perform proficiency testing with higher confidence.




- Test #
- Lot/Batch #
- Shelf-life Information
- **Truecal: Lot Specific Calibration**



Accurate Results – Impact of Truecal

Example: LCK303 High Range Ammonium



- Reduces variation
- Increases accuracy of results
- Higher Confidence

Compatible with DR3900 and DR6000

Focus on WW parameters: COD (LCK014 / 1414 / 114 / 314 / 514 / 614 / 714 / 914 / 1014) - Total N (LCK138 / 238 / 338 / 438) - Ammonium (LCK302 / 303 / 304 / 305 / 503 / 505) - Nitrate (LCK339 / 340 / 540) - Nitrite (LCK341 / 342) - Phosphate (LCK348 / 349 / 350)



Accurate Results – Always up-to-date with RFID

Reliable updates in a flash



Step 1

The DR 3900 reads the barcode on the cuvette and identifies that a data update is required for the test in question.



Step 2

It then indicates that the cuvette packaging should be placed in front of the photometer. After two seconds, an audio signal confirms that the data has been updated.



Step 3

Measurement is initiated automatically once the update is complete — using the correct factors to make sure that the measurement result is also correct.



Without RFID sensor the method can be updated via file download from www.hach.com

Accurate Results – Addista



Auswertformular • Evaluation form • Formulaire d'évaluation • Utvärderingsformulär
 Uitwaarderingsformulier • Modulo valutazione dati • Formulario de evaluación

	I	II	III	IV	V	VI
Küvetten-Test Cuvette-Test Test en Cuve Kyveti-test Kuvetten-test Cuvette-test Cubeta-Test	LCK1414	LCK313	LCK335	LCK311	LCK353	LCK341

Chargen-Nr.:
Batch No. of
Número do lot.
Provinummer:
Chargennummer:
n° lotte.
Número de lote.

ADDISTA: LCA 700

Chargen-Nr.:
Batch No. of
Número do lot.
Provinummer:
Chargennummer:
n° lotte.
Número de lote.

18176

Ringversuchlösung A:
Ringtest solution A:
Solution pour l'essai inter-laboratoires A:
Ringlösungslosning A:
Ringonderzoek-oplossing A:
Soluzione Ring-test A:
ES: Solución ringtest A

1 21,5 m

2

Ringversuchlösung B:
Ringtest solution B:
Solution pour l'essai inter-laboratoires B:
Ringlösungslosning B:
Ringonderzoek-oplossing B:
Soluzione Ring-test B:
ES: Solución ringtest B

1 27,5 m

2

Standardlösung:
Standard solution:
Solution standard:
Standardlösning:
Standaardoplossing:
Soluzione standard:
Solución estándar

1 30,7 m

2

HACH
Be Right™

Monsieur Rocton
ETS Mermier Lemarchand

LCK 313
Chrome
gamme de mesure 0,05 - 1,0 mg/l
Jour d'entree 13.12.2018
instrument de mesure LASA 100

Essais Inter-Laboratoires conformes ISO 572


Lösung A

Lösung B

Standard

	votre résultat		vaieur prévue
	[mg/l]		[mg/l]
Lösung A	0,187	0,186	0,19
Lösung B	0,445	0,441	0,45
Standard	0,499	0,502	0,5

Several parameters in the same standard solution



Addista-Standards for system control and free participation in HACH round robin test

Traceable Multi-Standards for accredited labs



Accurate Results – ISO Standards

Certificate of Conformity

We hereby declare that the
LANGE Cuvette Test LCI 400
is completely conform to
ISO 15705
„Determination of the chemical oxygen demand index (ST-COD) –
Small-scale sealed tube method“

Requirements of ISO 15705	Specifications of LCI400
Formulation: conc. (K ₂ Cr ₂ O ₇) = 4.596 g/L conc. (Ag ₂ SO ₄) = 9.375 g/L conc. (HgSO ₄) = 25.00 g/L conc. (H ₂ SO ₄) = 1435.2 g/L	conc. (K ₂ Cr ₂ O ₇) = 4.596 g/L conc. (Ag ₂ SO ₄) = 9.375 g/L conc. (HgSO ₄) = 25.00 g/L conc. (H ₂ SO ₄) = 1435.2 g/L
Sample Volume : min 2.0 mL	2.0 mL
Digestion : Temperature = 150 +/- 5 °C Time = 2 h +/- 10 min	Temperature = 148 °C Time = 2 h
Photometric Detection : Wavelength = 600 +/- 20 nm	Wavelength = 605 nm

Dr. Bruno Lange GmbH & Co. KG
QA Management
Düsseldorf, June 2003

Dr. Ralf Kloos

LANGE
FOR WATER QUALITY

DR. BRUNO LANGE
GMBH & CO. KG
Wissenschaftsstraße 11
D-40249 Düsseldorf
Tel. +49 (0)211 5288-0
Fax +49 (0)211 5288-143
kundenservice@brunolange.de
www.brunolange.de

COD ISO 15705:2002


INTERNATIONAL STANDARD **ISO 23697-1**

Water quality – Determination of
total nitrogen (TN_b) in water using
small-scale sealed tubes –
Part 1:
Dimethylphenol colour reaction

RECIPIENTS OF THIS DRAFT ARE INVITED TO SUBMIT WITH THEIR COMMENTS, NOTIFICATION OF ANY RELEVANT PATENT RIGHTS OF WHICH THEY ARE AWARE AND TO PROVIDE SUPPORTING DOCUMENTATION.

IN ADDITION TO THEIR EVALUATION AS BEING ACCEPTABLE FOR INDUSTRIAL, TECHNOLOGICAL, COMMERCIAL AND USER PURPOSES, DRAFT INTERNATIONAL STANDARDS MAY ON OCCASION HAVE TO BE CONSIDERED IN THE LIGHT OF THEIR POTENTIAL TO BECOME STANDARDS TO WHICH REFERENCE MAY BE MADE IN NATIONAL REGULATIONS.

Reference number
ISO/FDIS 23697-1:2022(E)



© ISO 2022

Ammonium: ISO 23695

TN: ISO 23697-1


Nitrate: ISO 23696-1



Environmental Protection - Recycling

Waste Handling in Düsseldorf

- ~75% of test components returned from our customers to the environment center will be re-used or recycled.
- Certified disposal company & registered as an official recycling center
- Commitment to protect the environment (ISO14001 certified) and a safe working place (ISO45001 certified)
- New technologies which enhance the process are actively implemented if an environmental benefit is found.



Deutscher Nachhaltigkeitspreis



COD



Be Right™

COD - Chemical Oxygen Demand

What is COD?

- Chemical Oxygen Demand
- Sum parameter, determined from a homogenized sample
- Organic compounds are oxidized
- Method derived parameter: Defined conditions $K_2Cr_2O_7$, 148 °C and 2 h digestion silver-catalysed in strong sulfuric acid)
- Determination by photometry (e.g. ISO 15705) or titration (e.g. ISO 6060) with $(NH_4)_2Fe(SO_4)_2$



ISO 6060 TOC (macro size, titration based)



ISO 15705 ST COD with 95% less reagents (micro size, sealed tube COD, photometric based)


COD – Sample Preparation

- **Homogenization of sample is mandatory**, according standards (e.g. DIN 38402-30)
- Particles, fibers or flakes (precipitation) are part **sum parameter**
- **Representative sample**
- **Dispersing tool** needed for homogenization
- Specification: **20,000 rpm for 60 sec.**
- Longer time frames **can influence** the sample
- Final **sample taking** for analysis purposes with **magnetic stirrer**



COD – Procedure

Precise COD measurement in six steps



- Shake cuvette before adding sample
- Add sample into “cloud”
- Close cuvette, start digestion
- Shake warm cuvette again after digestion and cool down
- Make sure the precipitation has settled completely before reading and do not swirl cuvettes just before reading
- Start reading on photometer

Note

- The color is stable (Cr^{3+}) in the cuvette
- COD tests readings can be done next morning



COD – Chloride Interference



- Chloride in **high concentrations** lead to false positive results, check **interference – table** in each test manual!
- **Check chloride** content in sample prior to COD analysis (e.g. LCK 311, Chloride test strips)
- If the **chloride concentration is above the interference limit** of e.g. **1500 mg/l**, sample must be diluted with **COD-free water**

or better: apply COD tests with higher Chloride tolerance:

- LCK 1014 COD (100-2000 mg/l), Cl⁻ Tolerance 4,000mg/l
- LCK1714 COD salt water (70 – 250 mg/L), Cl⁻ Tolerance 20,000 mg/l
- LCK 1814 COD salt water (7 - 70 mg/L), Cl⁻ Tolerance 20,000 mg/l
- LCK 1914 COD salt water (250 - 1000 mg/L), Cl⁻ Tolerance 20,000 mg/l

COD – Choosing the correct range

- Samples with COD values higher than measuring range can be diluted, but usually **not preferred** for COD measurements
- **Different ratio** between dichromate and COD concentration in cuvette, different **oxidation performances**
- Highest accuracy and precision **in the middle of the range**
- **Different test ranges** available to avoid **inaccuracies**



5 – 60 mg/l



15 – 150 mg/l



50 – 300 mg/l



150 – 1000 mg/l



100 – 600 mg/l



100 – 2000 mg/l




1000 – 10000 mg/l



5000 – 60000 mg/l



- Example: COD is 145 mg/l (LCK 314: 15 – 150 mg/L)
- Better to choose LCK 614 (50 – 300 mg/l)



COD – High Temperature Application

HT COD application:

- Companies that **need fast COD results** (examples):
 - unload waste truck
 - fast results for plant control
- HT 200S allows COD digestion at **170°C in 15 minutes**
- **HT application:** If sample is **unknown** start first with **comparative trials** (HT \
- **Municipal wastewater** shows in general **good matches** to standard (ISO 15705, 148°C, 2h)
- For correct documentation the results can be **marked** as **HT-COD**



COD – Standard Methods



Requirements of ISO 15705 Formulation (after volume contraction):	Specifications of LCI 400 (after volume contraction)
conc. (K ₂ Cr ₂ O ₇) = 4.845 g/L	conc. (K ₂ Cr ₂ O ₇) = 4.845 g/L
conc. (Ag ₂ SO ₄) = 9.881 g/L	conc. (Ag ₂ SO ₄) = 9.881 g/L
conc. (HgSO ₄) = 26.35 g/L	conc. (HgSO ₄) = 26.35 g/L
conc. (H ₂ SO ₄) = 1512.7 g/L	conc. (H ₂ SO ₄) = 1512.7 g/L
Sample Volume: min. 2.0 mL	2.0 mL
Digestion: Temperature = 150 °C +/- 5 °C Time = 2 h +/- 10 min	Temperature = 148 °C Time = 2 h
Photometric Detection: Wavelength = 600 nm +/- 20 nm	Wavelength = 605 nm

HACH LANGE GmbH
LSBU Application
Düsseldorf, February 2011

R. Kloos
Dr. Ralf Kloos





HACH LANGE
GMBH
Wilschstraße 11
D-40249 Düsseldorf
Tel: +49 (0)211 2389-0
Fax: +49 (0)211 2389-143
kundenservice@hach-lange.de
www.hach-lange.de

Ministerieel besluit van 22 januari 2021 — Belgisch Staatsblad van 03 februari 2021

Compendium voor de monsterneming, meting en analyse van water



COD Saltwater Application LCK1714, 1814 1914 – Flemish Standard WAC/III/D/020

- New Flemish (Belgium) COD Standard: WAC/III/D/020 (1/2021) Vito (independent Flemish research organization)
- WAC/III/D/020 is mandatory to use for COD analytics in saltwater, launched Jan 2021
- Macro size reflux method (e.g. ISO 6060) not used for such samples any longer
- Only micro size Cuvette Tests
- Together with Vito, - Hach introduced a greener analysis method for measuring COD in wastewater (s. "Vision" newsletter, page 14-15):

Bepaling van het chemisch zuurstofverbruik (CZV)



Versie oktober 2020

WAC/III/D/020



COD – Regulatory Shift to TOC?

REACH

- **REACH directive:** 2012 ECHA (EU Commission) started activities to ban Cr⁶⁺, SVHCs (“Substance of very high concerns”)
- **Conflict:** COD as **wastewater** parameter, COD **requested** by UWWTD
- **Authorization** needed to **keep working** with Cr⁶⁺ and COD test **Sunset Date (09/2017)**
- ECHA **exempted** SVHC’s (Cr⁶⁺), if used for **scientific research & development** (monitoring of water) and in **amounts < 1t/a**

Current state:

- ~ **2016:** DG Environment (EU Commission) **seeking for alternatives** to COD (go green strategy)
- Several drafts to replace COD by TOC short term
- Plant **specific conversion factors** suggested, for each WWTP (Waste Water Treatment Plant)
- **General replacement not accepted**, but **total revision** of UWWTD started
- Some Member States, e.g. NL start investigations to **replace COD by TOC**

Future state:

Revised EU UWWT Directive October 2022

- TOC official monitoring parameter, member states can choose between TOC or COD
- Lower discharge limits for TN and TP
- Removal and Monitoring of Micropollutants (e.g. painkiller, antibiotics, anti-corrosion agents...)
- Wastewater surveillance for multiple viruses (incl SARS-CoV- II)
- NL moves away from COD to TOC in 2025 for all sewage charge testing





TOC




Be Right™

What is TOC?



KHP Potassium hydrogen phthalate



Calcium Carbonate

TOC Detection methods:
1) Direct Method
2) Difference Method

TOC Direct Method




LCK385 - LCK387: Persulfate Oxidation after IC removal

IC removal

- Sample + diluted sulfuric acid -> vibration and fan in TOC-X5 shaker
- 5 min vibration sufficient to remove 250mg/L TIC

Mechanism

- persulfate digestion 100° C / 120 min
- Oxidation of organic carbon
- CO₂ passes gas permeable membrane
- Color change of pH indicator
- Ranges: 3-30mg/L; 30-300mg/L;300-3000 mg/L



TOC Difference Method

LCK380-LCK381: difference method

- For samples containing VOC or low TIC
- 2 cuvettes required: 1 for TC and 1 for TIC
- Photometer calculates
 $TOC = TC - TIC$





Phosphate



Be Right™

Phosphate Overview



Phosphate – Regulatory Outlook

Water Framework Directive (WFD)

- WFD's goal is to **increase water quality** by 2027
- Achieving **“good status”** for rivers, lakes, coastal waters and groundwater.
- Wastewater Treatment Plants (WWTP's) **“point-sources”**

Outlook:

- Ongoing **lowering of Phosphorus concentrations** in surface waters **until 2027 to 0.05 mg/L P**
- **Challenging treatment** of existing WWTP's
- **Precipitation methods** of WWTP's can reach their limits
- **New treatment techniques & new monitoring** systems to come



Phosphate – Hach solutions addressing WFD goals

Laboratory based measurement

Ultra-low range Phosphate

Range: **0.01-0.5mg/l PO₄-P** (ortho)

50mm cuvette with Dosiscap

DR3900/DR6000

Phosphor molybdenum blue-method

Method according ISO 6878-1-1986



Continuously measurement solution

Phosphax sc LR (Low Range)

Range: **0.015-2 mg/L PO₄-P** (ortho)

Yellow method

Aeration tank and WWTP effluent


Surface water monitoring



Claros Process Management

Claros RTC (Real Time Control)

RTC-P module for Phosphorous removal



Phosphate – Different options

- Phosphate tests to measure different phosphate components
- Different sample **preparation & digestion** of sample
- Evaluation as **ortho PO₄ –P** or **TP** (total Phosphorous)
- Legal reporting: **WW** as **PO₄-P**, **DW** as **PO₄**



Target	Sample	Oxidation	Digestion
Total Phosphate	Homogenization	Yes	Yes
Ortho-Phosphate	Filtration	No	No
Polyphosphate (Ortho-phosphate is co-determined)	Filtration	No	Yes
Phosphonates (O-Phosphate & Polyphosphates are co-determined)	Filtration	Yes	Yes



Phosphate – Procedure

- Analysis done out of **homogenized** sample:

- Shake (up to 5 l)
- Magnetic stirrer
- Disperser
- Ultrasonic bath



- Choose correct **measuring range** of test
- Otherwise, dilute sample

- Oxidation in Reactor:**

- HT200S: 15 mins
- LT200: 30 mins, 120°C
- LT200: 60 minutes at 100° C




➔ Higher Temperature to speed up digestion time


Procedure total phosphorus




1. Carefully remove the foil from the screwed-on DosiCap Zip.




2. Unscrew the DosiCap Zip.




3. Carefully pipet 2.0 mL sample.




4. Immediately screw the DosiCap Zip back on tight; fluting at the top.




5. Shake vigorously.




6. Heat in the thermostat. HT 200 S; in the standard program HT for 15 minutes. Thermostat: for 60 minutes at 100° C (212° F) or for 30 minutes at 120° C (248° F).




7. Allow to cool to room temperature. NOTE: Check if the cap is still tight after cooling.



8. Shake vigorously.




9. Unscrew the DosiCap Zip.




10. Pipet into the cooled cuvette: 0.2 mL Reagent B. Close Reagent B immediately after use.




11. Screw a grey DosiCap C on the cuvette.



12. Invert a few times until the freeze-dried contents are completely dissolved.



13. After 10 minutes, invert a few more times, thoroughly clean the outside of the cuvette and evaluate.




14. Insert the cuvette into the cell holder. DR1900: Go to LCK/TNTplus methods. Select the test, push READ.

Nitrogen Parameters




Be Right™


Nitrogen – Definition Nitrogen Parameters




Organic Nitrogen (N_{Org})



Urea:
Important role in biological processes, e.g. metabolism of proteins



Protein (protein):
Biological macromolecule, formed from amino acids by peptide bonds



Organic acids (amino acids):
e.g. Tyrosine, non-essential proteinogenic α-Amino acid, casein (cheese)

Total Kjeldahl-Nitrogen TKN

- Total Kjeldahl-Nitrogen:**
- Sum of Ammonium + Organic Nitrogen
 - Analyzed by "Kjeldahl-Method"



Ammonium

- ISO standard (ISO 23695)
- Ammonium ions are not bounded/soluble, direct readings
- Filtering if turbidity is present
- Highest precision/accuracy in the middle of MR
- Multiple ranges enable most accurate testing without diluting, e.g. Ammonium

Quality certificate
 Technical data for Validation
 of LCK305 (1.3-15 mg/l
 Ammonium)

Quality certificate

Technical data for cuvette test LCK305
 (Results as NH₄)

Sensitivity	0.140 Abs./(mg/l)
Ordinate Intercept	0.092 Abs.
Residual standard deviation	0.0019 Abs.
Method variation coefficient	1.63 %
Method standard deviation	0.13 mg/l
Confidence Interval (95%)	± 0.33 mg/l


The technical data for cuvette test LCK305 were determined in conformity with ISO 8466-1 and DIN 38402 A51 „Calibration of analysis methods“.

The series of the smallest and largest calibration standards exhibit normal distribution and are outlier- and trend-free. The calibration gives a linear function.

Technical data in conformity with DIN 32645

Detection limit	0.033 mg/l
Quantitation limit	0.099 mg/l

The detection and the quantitation limits were determined in conformity with DIN 32645.




Result	Confidence Interval (95%)
3.0 mg/l	± 0.35 mg/l
6.0 mg/l	± 0.33 mg/l
9.0 mg/l	± 0.33 mg/l
12.0 mg/l	± 0.34 mg/l
15.0 mg/l	± 0.36 mg/l

HACH LANGE GmbH
 Quality Management


R. Kloos
 Dr. Ralf Kloos



HACH LANGE GmbH
 Wilstättstraße 11
 D-40549 Düsseldorf
 Tel. +49 (0)211-52 88-320
 Fax +49 (0)211-52 88-143
 Kundenservice@hach-lange.de
 www.hach-lange.de



Ammonium LCK ranges in mg/L NH₄-N




Ammonium – Tipps & Tricks

- Interferences from water sample are stated in each working manual
- If colored sample or turbidity is present, a separate determination of a **sample-specific-blank** according Hach's Application A 130 might be required

Test	Cuvette	Sample	Dist. water	Cap	Reagent A	Reagent B
LCK 302 Ammonium	Analysis cuvette	0.2 mL	—	Original cap	—	—
LCK 303 Ammonium	Analysis cuvette	0.2 mL	—	Original cap	—	—
LCK 304 Ammonium	Analysis cuvette	5.0 mL	—	Original cap	—	—
LCK 305 Ammonium	Analysis cuvette	0.5 mL	—	Original cap	—	—
LCK 307 Boron	Empty cuvette*	2.5 mL	—	Red plug	1.0 mL	—
LCK 310 Chlorine	Empty cuvette*	2.0 mL	—	Red plug	—	—

- The chemical reaction is temperature-dependent (room temperature)
- If temperature is too warm (e.g. non-air conditioned laboratory in mid-summer): **Application „Temperature Compensation“** is recommended (only DR3900/6000)
- The chemical reaction is very time-dependent (15 min) and must be strictly observed



Sample Specific Blank Value
Detecting and compensating for turbidity and colouration

During a photometric analysis, in addition to the compounds listed in the interferences table turbidity and colouration can lead to false results. These errors come either from the sample itself or arise from reactions with the reagents. The influence of the turbidity and/or colouration can be determined by establishing the sample specific blank value.

Contents

- What is a sample specific blank value?
- Working procedure / measurement / evaluation
- Check list of sample specific blank value
- Special evaluation using blank value correction

APPLICATION NOTE: Temperature Compensation for the LCK cuvette tests Ammonium (LCK303, LCK304 and LCK305)

Temperature compensation for the LCK cuvette tests Ammonium (LCK303, LCK304 and LCK305)

Dr. Ralf Kloos, Product Application Manager

Introduction

The LCK cuvette tests for Ammonium are based on a reaction to form indophenol. Ammonium ions react at pH 12.6 with hypochlorite ions and salicylate ions in the presence of sodium nitroprusside as a catalyst to form indophenol. The amount of color formed is directly proportional to the amount of Ammonium Nitrogen present in the sample. The color development is measured at 694 nm and the results are expressed as mg/L NH₄⁺-N or NH₄⁺.

The LCK cuvette tests are calibrated by a ten point-calibration according to ISO 8466-1 at a defined temperature of 20°C. The reaction described above is temperature-sensitive, so the recommended sample and reagent temperature is 20°C. Incorrect results may be obtained if the test is not performed at the recommended temperature.

Hach has performed the tests at defined temperatures of 15°C-30°C and calculated the bias to the results obtained at 20°C. The results of this evaluation are shown in the following Application Note.

Material

- LCK cuvette tests for Ammonium - LCK303, LCK304, and LCK305
- Hach spectrophotometers DR3900 and DR6000
- JULABO water bath assembly to achieve defined temperatures from 15°C – 30°C
- Hach Ammonia Standard Solution - 1000 mg/L as NH₄⁺-N Cat. 23541-53
- Hach calibrated pipettes (0.2-1.0 mL and 1.0-5.0 mL) and pipette tips

Nitrate – Tipps & Tricks

- ISO standard (ISO 23696)
- Reliable results in 5 steps (example LCK340):


Procedure



1. Carefully pipet 0.2 mL of sample.




2. Carefully pipet 1.0 mL of solution A.



3. Close the cuvette and invert a few times until no more streaks can be seen.



4. After 15 minutes, thoroughly clean the outside of the cuvette and evaluate.



5. Insert the cuvette into the cell holder.

- Perform analysis in a timely manner after sampling or preserve the sample
- Not more than 3 hours should elapse between sampling and analysis.
- High loads with **oxidizable organic substances (COD > 200 mg/l)** lead to **additional coloration** of the reagent (false positive results)
- Chloride concentrations > 500 mg/L in a samples interfere (false negative results)
- If dilution is no option due to too low NO₃ concentration, Chloride removal using Silver Oxide (LCW 925, up to 20 g/l Cl⁻) is recommended



Total N – Tipps & Tricks

Principle:

- ISO standard (ISO 23697)
- **Inorganic and organically bound Nitrogen** is oxidized to Nitrate
- Thermal digestion in **reaction tubes** with Peroxydisulfate
 - 30 minutes at 120 °C (LT200)
 - 15 minutes with HT200s
- **Photometric determination** and calculation to TNb

Interferences:

- **Sample homogenization** is mandatory
- Dependent on TN range, specific concentrations of chloride or high COD may interfere with the digestion and result in low findings (see test procedures!)
- Reaction tubes (20 mm) recommended to use **max 7 times** (carry over effects, **Nitrogen-free** water needed)
- Slight **turbidities present** do not interfere
- During the reaction, a slight pink coloration may form. The coloration does not affect the evaluation.



1-16 mg/L TNb



5-40 mg/L TNb



20-100 mg/L TNb



New range available! LCK438 100-250 mg/L TNb

Thank you!

