

TOC Too Hard? No, Too Easy With Hach!

Use Hach's latest colorimetric test for TOC: the TNTplus TOC reagent and TOC-X5 shaker used alongside Hach's easy-to-use spectrophotometers.

**30%
OFF**



Save up to 60%
of your time for TOC
testing!

Advantages Versus Other Colorimetric TOC Measurements

- Fast and easy 5-minute processing time to remove total inorganic carbon (TIC)
- Auto detection and recognition with barcode on vial
- No blank needed so there is less waste
- No handling and breaking glass indicator vials

Advantages Versus Dedicated TOC Lab Analysers

- No need for a large capital outlay since it uses a spectrophotometer
- Easy to use: follow the instructions on the box

**30%
OFF**



Buy the TOC-X5 & TNT810-AU
and **Get 30% Off***



Interested in buying a new Hach
DR3900 or DR6000 Spectrophotometer?
Speak to your Sales Representative!

Terms & Conditions

- Get 30% discount off the list price on TOC-X5 Shaker and TNTplus TOC reagent from 1st October 2024 to 20th December 2024
- TOC-X5 Shaker only compatible with Hach Spectrophotometers and TOC TNTplus reagents.
- Cannot be combined with other Discounts or Offers.
- Not intended for Distributors.
- Not valid for Orders in the past.
- Valid for AU and NZ customers only.



www.au.hach.com

Contact us for more info: (Australia)
1300 887 735 | sales@hachpacific.com



Be Right™